PRODUCT

VISUCORE® 500
by ZEISS

INSTRUCTIONS FOR USE

PRODUCER

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
OPERATING CONDITIONS	
Relative humidity:	10-90 %
Operating temperature:	+10°C - +40°C
Pressure:	500-1060 hPa

ACCOMPANYING DOCUMENTS	
For full parts list and specification, please refer to the technical specification document	

OPERATION

Activation

When supplied with electricity, the Table is active.
Visually inspect cables and electrical components for damage. Never use the product if you spot any damaged electrical components.

Height adjustment

Teburu XR ONE SC table is equipped in two hand switches and one foot switch to adjust the height of the table. Hand and foot set incorporates two buttons with triangular symbols corresponding to 'Up' and 'Down'. Press and hold button until the tabletop attains a desired height.
SAFETY PRECAUTIONS: Please avoid any hazards to the patient during the height adjustment.

Deactivation

When disconnected from the electricity, the Table is inactive.

Wheels locking

Teburu XR ONE table has a central locking system of wheels. You can lock all four wheels simultaneously with foot pedals located in both the cranial and caudal end of the table.

FAQ

Cleaning and maintenance

Use only soft cleaning agents and high quality cloth. To avoid bruising and smudging, polish the recently cleaned component with a soft, dry cloth. Please avoid soaking electricity components (cords, inlets etc.)

Disposal

Please see that Teburu fixation table contain parts that cannot be disposed in domestic waste. Please contact a waste management expert should you want to dispose of these.

WARNING AND SAFETY NOTICES

Intended and Unintended applications

Teburu XR ONE SC is a mobile fixation table which gives excellent possibilities to index all different kinds of immobilization devices on top of it. Other applications than specified may cause hazards and/or damage.

Liability of the Manufacturer

CX Engineering Sp. z o.o. considers itself responsible for issues on safety, reliability and performance only when:

- Any modifications and repairs are conducted by authorized persons AND.
- The product is used according to the manual.

Responsibility of the User

- Ensuring that electric supply is correct (see table on page 1).
- Cleaning and maintenance.
- Operation according to the manual.
- Storage in safe conditions.

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Teburu table may be individually built. To recognize the different versions there is a code in a format that indicates the type of the table TEBURU system, when any letter from the table, attached to number from 0 to 99999 may stand for: „L”, „R”, „L”, „R”, „L”, „R”, „L”, „R”. Additionally, any character from this code may be removed if table configuration does not require it.

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Scope: Design and development, production, distribution and service of medical furniture

The Certification Body of TÜV Rheinland LGA Products GmbH certifies that the organization has established and applies a quality management system for medical devices. Proof has been furnished that the requirements specified in the documented standard are fulfilled. The quality management system is subject to yearly surveillance.

Report No.: 84948524-30
Effective date: 2021-01-19
Expiry date: 2023-10-13
Issue date: 2021-01-19

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EXPECTED POSITIONS

When using the table, Users are expected to be standing within safe distance when the height is adjusted. If no adjustment of the table setting is required, the wheels must be locked. Only then the users may use the tabletop to support their elbows or hand or palms.

EXPECTED POSITIONS

Please make sure that:

- the table is unplugged from electricity,
- no cords or plugs are hanging loose,
- the column is refracted to the shortest possible length,
- the equipment installed on the table does not exceed the contour of the table (e.g. edges of the table top, edges of the base plate),
- when pushing over thresholds the table is positioned sideways (i.e. longer edges of the tabletop are perpendicular to the obstacle),
- when pushing or dragging the table the force is not to be applied to the handset, cords, plugs,
- no body, clothing or equipment parts are to be positioned between the contour of the table and obstacles such as walls, door frames and other, when the distance between those shortens to values that may result in harm or hazard due to squeezing,
- feet are not to be run over by table wheels.